

No matter where you are testing, Benson Medical Instruments provides seamless solutions for industrial hearing testing, earplug fit testing, and spirometry testing for hearing conservation and pulmonary function program management.

The full line of stand-alone and computer-controlled audiometers, fit testers, and spirometer delivers fast audiograms and spiograms, fully integrated databases with network connectivity, tools to handle difficult-to-test subjects, increased reliability, and unsurpassed ease in handling data. Benson Medical products are designed to meet OSHA/MSHA/DOT, ATS/ERS and NIOSH/CWHSP standards.

Audiometers



CCA-200mini Audiometer – Your most powerful option - Ideal for testing of populations over 200 or managing data with Plus Package for multiple companies.

CCA-100mini Audiometer – Ideal for testing 1 - 2 stations.

Database Integration



Spirometry



CCS-200 Computer Controlled Spirometer
Produces accurate and repeatable spiograms for use in NIOSH-compliant programs. Exceptional hygiene provided by single-use "U-tube" mouthpiece. Integrates with the CCA-200mini, CCA-100mini audiometers, and CCF-200 Fit Test.

Fit Testing



CCF-200 Fit Test – Automatically measures the sound attenuation that a subject receives while wearing earplugs. Provides single or multiple station testing. Works with a wide variety of manufacturers' earplugs. Database integrates with the CCA-200mini and CCA-100mini audiometers, and CCS-200 spirometer.

NEXT Audiometer



The ultimate stand-alone audiometer – touchscreen operation with powerful analysis for OSHA compliance.

BAS-200 & BAS-200slm



Bio-acoustic simulators that perform dependable OSHA-required daily calibration check of Benson audiometers. The BAS-200slm with sound level meter allows for OSHA and ANSI testing without a sound room.

Solo™ Software



Allows you to transfer audiograms from your third-party audiometer for analysis and reporting. Calculates STS's and recordable shifts, and tracks medical referrals for single or multiple companies.

Benson Cloud Services



Integrates with the CCA-200mini and CCA-100mini audiometers to deliver easy-to-use audiogram management. Secures audiogram data in a HIPAA-compliant cloud-based data management system. Access your data via browser with a password. Always get automatic data backup and the most current features.

BENSON MEDICAL INSTRUMENTS

310 Fourth Avenue South
Suite 5000

Minneapolis, MN 55415

Learn more at

www.bensonmedical.com or

by contacting

Benson Medical Sales:

sales@bensonmedical.com

612-827-2222 Ext. 3